



**Wellmien Healthcare Tech. (Suzhou) Co., Ltd.**

3603, Global 188, No.88 Suhui Rd, Suzhou, Industrial Park,  
Jiangsu, China, 215021

Tel : +86 512 65953508

Web: [www.wellmien.com](http://www.wellmien.com)

[www.wellmien.com.cn](http://www.wellmien.com.cn)

## Technical Data Sheet


	<b>Item</b>	Disposable PE Apron
	<b>Material</b>	polyethylene
	<b>Weight</b>	35g+-2g
	<b>Color</b>	Blue
	<b>Style</b>	Long sleeve with thumb hook design, 2 ties at waist, perforated back for easy removal
	<b>Size</b>	See below
	<b>Packing</b>	1pc/bag 10bags/inner bag, 200pcs/carton



### Container Loading

Carton Size (cm)	Gorss Weight	Pcs carton	Ctn/ 20FT	Pcs/ 20FT	Ctn/ 40HQ	Pcs/ 40HQ
41*27*35	4KG	200	725	145,000	1400	280,000

	Spec (cm)	Tolerance (cm)		Spec (cm)	Tolerance (cm)
A	190	+/-2.5	H	53	+/-2.5
B	115	+/-1.3	I	13	+/-1.3
C	12.5	+/-1.3	J	75	+/-2.5
D	82	+/-1.3	K	35	+/-2.5
E	32	+/-1.0	L	75	+/-2.5
F	118	+/-1.3			
G	8.5	+/-1.3			





**Wellmien Healthcare Tech. (Suzhou) Co., Ltd.**


3603, Global 188, No.88 Suhui Rd, Suzhou, Industrial Park,  
Jiangsu, China, 215021

Tel : +86 512 65953508

Web: [www.wellmien.com](http://www.wellmien.com)

[www.wellmien.com.cn](http://www.wellmien.com.cn)

## KEY FEATURES



## INTENDED USAGE:

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.
- Food processing for personal protection

## STORAGE And TRANSPORT:

- Keep in ventilated and dry room, prevent heavy weight, direct sunlight, rain and snow immersion.
- Shelf Life: 5 years.